

## COOLEY SUCTION TUBE

Ø 7.0mm; 30.0cm (12")

**SKU:** 085-024-03

**Categories:** [Open Surgery](#), [Puncture / Special Needles](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Suction Tubes</a>
<b>KLS Martin</b>	18-565-30
<b>AUTHORS</b>	<a href="#">COOLEY</a>
<b>Medicon</b>	09.86.15

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Suction Tubes</a>
<b>KLS Martin</b>	18-565-30
<b>AUTHORS</b>	<a href="#">COOLEY</a>
<b>Medicon</b>	09.86.15