

## COOLEY VASCULAR DILATORS

flexible; Ø 3.0mm

13.0cm (5.2")

**SKU:** 141-316-30

**Categories:** [Open Surgery](#), [Probes](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vascular Dilators</a>
<b>KLS Martin</b>	24-490-30
<b>AUTHORS</b>	<a href="#">COOLEY</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vascular Dilators</a>
<b>KLS Martin</b>	24-490-30
<b>AUTHORS</b>	<a href="#">COOLEY</a>