

COOLEY VASCULAR DILATORS

flexible; Ø 3.0mm

13.0cm (5.2")

SKU: 141-316-30

Categories: [Open Surgery](#), [Probes](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|-----------------------------------|
| UNIT | pcs. |
| TYPE | Vascular Dilators |
| KLS Martin | 24-490-30 |
| AUTHORS | COOLEY |

ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|-----------------------------------|
| UNIT | pcs. |
| TYPE | Vascular Dilators |
| KLS Martin | 24-490-30 |
| AUTHORS | COOLEY |