

## COTTLE RETRACTOR

2-pronged; sharp

14.5cm (5.8")


**SKU:** 112-308-00

**Categories:** [Open Surgery](#), [Rhinology](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Delicate Hooks</a>
<b>KLS Martin</b>	37-708-14
<b>AUTHORS</b>	<a href="#">COTTLE</a>
<b>SPECIFICATIONS</b>	<a href="#">sharp</a> , <a href="#">toothed</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Delicate Hooks</a>
<b>KLS Martin</b>	37-708-14
<b>AUTHORS</b>	<a href="#">COTTLE</a>
<b>SPECIFICATIONS</b>	<a href="#">sharp</a> , <a href="#">toothed</a>