

COUPLAND TUBE ADAPTER

for COUPLAND inserts; 16.5cm

SKU: 085-010-17

Categories: [Open Surgery](#), [Puncture / Special Needles](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Suction Tubes
KLS Martin	18-575-24
AUTHORS	COUPLAND

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Suction Tubes
KLS Martin	18-575-24
AUTHORS	COUPLAND