

## COUPLAND TUBE ADAPTER

for COUPLAND inserts; 16.5cm

**SKU:** 085-010-17

**Categories:** [Open Surgery](#), [Puncture / Special Needles](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Suction Tubes</a>
<b>KLS Martin</b>	18-575-24
<b>AUTHORS</b>	<a href="#">COUPLAND</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Suction Tubes</a>
<b>KLS Martin</b>	18-575-24
<b>AUTHORS</b>	<a href="#">COUPLAND</a>