

CRILE NEEDLE HOLDER

straight; serrated; without catch

15.0cm (6")

SKU: 050-338-15

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM503
KLS Martin	20-500-15
AUTHORS	CRILE
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM503
KLS Martin	20-500-15
AUTHORS	CRILE

SPECIFICATIONS

[serrated](#), [straight](#)