

## CUSHING RETRACTOR

fenestrated

20.5cm (8.2")


**SKU:** 062-250-20

**Categories:** [Open Surgery](#), [Retractors](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Delicate Hooks</a>
<b>KLS Martin</b>	15-170-20
<b>AUTHORS</b>	<a href="#">CUSHING</a>
<b>Medicon</b>	57.52.05



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Delicate Hooks</a>
<b>KLS Martin</b>	15-170-20
<b>AUTHORS</b>	<a href="#">CUSHING</a>
<b>Medicon</b>	57.52.05