

## CW-CASPAR DISC RONGEUR

straight; 3 x 12mm; toothed

18cm (7")



**SKU:** 240-511-03-CW

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Laminectomy Rongeurs</a>
<b>AUTHORS</b>	<a href="#">CASPAR</a> , <a href="#">CLEANWAVE</a>
<b>SPECIFICATIONS</b>	<a href="#">straight</a> , <a href="#">toothed</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Laminectomy Rongeurs</a>
<b>AUTHORS</b>	<a href="#">CASPAR</a> , <a href="#">CLEANWAVE</a>
<b>SPECIFICATIONS</b>	<a href="#">straight</a> , <a href="#">toothed</a>