

CW-CASPAR DISC RONGEUR

straight; 4 x 14mm; toothed

18cm (7")



SKU: 240-511-04-CW

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Rongeurs
AUTHORS	CASPAR , CLEANWAVE
SPECIFICATIONS	straight , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Rongeurs
AUTHORS	CASPAR , CLEANWAVE
SPECIFICATIONS	straight , toothed