

CW-CASPAR DISC RONGEUR

angled up; 2 x 12mm; toothed

18.0cm (7.4")

SKU: 240-515-02-CW

Categories: [Neurosurgery](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Rongeurs
AUTHORS	CASPAR , CLEANWAVE
SPECIFICATIONS	angled up , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Rongeurs
AUTHORS	CASPAR , CLEANWAVE
SPECIFICATIONS	angled up , toothed