

## CW-CASPAR DISC RONGEUR

angled up; 2 x 12mm; toothed

18.0cm (7.4")



**SKU:** 240-515-02-CW

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Laminectomy Rongeurs</a>
<b>AUTHORS</b>	<a href="#">CASPAR</a> , <a href="#">CLEANWAVE</a>
<b>SPECIFICATIONS</b>	<a href="#">angled up</a> , <a href="#">toothed</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Laminectomy Rongeurs</a>
<b>AUTHORS</b>	<a href="#">CASPAR</a> , <a href="#">CLEANWAVE</a>
<b>SPECIFICATIONS</b>	<a href="#">angled up</a> , <a href="#">toothed</a>