

CW-CASPAR DISC RONGEUR

angled up; 3 x 12mm; toothed

18cm (7")



SKU: 240-515-03-CW

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Rongeurs
AUTHORS	CASPAR , CLEANWAVE
SPECIFICATIONS	angled up , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Rongeurs
AUTHORS	CASPAR , CLEANWAVE
SPECIFICATIONS	angled up , toothed