

CW-CASPAR DISC RONGEUR

angled down; 2 x 12mm; toothed

18.5cm (7.4")

SKU: 240-518-02-CW

Categories: [Neurosurgery](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Laminectomy Rongeurs |
| AUTHORS | CASPAR , CLEANWAVE |
| SPECIFICATIONS | angled down , toothed |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Laminectomy Rongeurs |
| AUTHORS | CASPAR , CLEANWAVE |
| SPECIFICATIONS | angled down , toothed |