

## CW-CASPAR DISC RONGEUR

angled down; 3 x 12mm; toothed

18.5cm (7.4")

**SKU:** 240-518-03-CW

**Categories:** [Neurosurgery](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

|                       |   |
|-----------------------|---|
| <b>UNIT</b>           | <a href="#">pcs.</a>                                  |
| <b>TYPE</b>           | <a href="#">Laminectomy Rongeurs</a>                  |
| <b>AUTHORS</b>        | <a href="#">CASPAR</a> , <a href="#">CLEANWAVE</a>    |
| <b>SPECIFICATIONS</b> | <a href="#">angled down</a> , <a href="#">toothed</a> |



### ZUSÄTZLICHE INFORMATIONEN

|                       |   |
|-----------------------|---|
| <b>UNIT</b>           | <a href="#">pcs.</a>                                  |
| <b>TYPE</b>           | <a href="#">Laminectomy Rongeurs</a>                  |
| <b>AUTHORS</b>        | <a href="#">CASPAR</a> , <a href="#">CLEANWAVE</a>    |
| <b>SPECIFICATIONS</b> | <a href="#">angled down</a> , <a href="#">toothed</a> |