

CW-CASPAR DISC RONGEUR

angled down; 3 x 12mm; toothed

18.5cm (7.4")



SKU: 240-518-03-CW

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Rongeurs
AUTHORS	CASPAR , CLEANWAVE
SPECIFICATIONS	angled down , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Rongeurs
AUTHORS	CASPAR , CLEANWAVE
SPECIFICATIONS	angled down , toothed