

CYSTOSCOPE HD+

Ø 4.0mm; 70°; 300 mm

autoclavable

SKU: 601-104-70HD

Categories: [Endoscopy](#), [Equipment](#), [Gynecology](#), [Urology](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
RANGE [Telescopes](#)
SCOPE [Gynecology](#), [Urology](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
RANGE [Telescopes](#)
SCOPE [Gynecology](#), [Urology](#)