

## CYSTOSCOPE HD+2

Ø 4.0mm; 12°; 300 mm

autoclavable

**SKU:** 601-104-12HD2

**Categories:** [Endoscopy](#), [Equipment](#), [Gynecology](#), [Urology](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Telescopes</a>
<b>SCOPE</b>	<a href="#">Gynecology</a> , <a href="#">Urology</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Telescopes</a>
<b>SCOPE</b>	<a href="#">Gynecology</a> , <a href="#">Urology</a>