

## DANDY HAEMOSTATIC FORCEPS

curved; serrated

14.0cm (5.6")



**SKU:** 040-345-14

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a> , <a href="#">Dissecting and Ligature Forceps</a>
<b>Aesculap</b>	BH139
<b>KLS Martin</b>	13-355-14
<b>AUTHORS</b>	<a href="#">DANDY</a>
<b>Medicon</b>	16.11.84
<b>SPECIFICATIONS</b>	<a href="#">curved</a> , <a href="#">serrated</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Artery Forceps</a> , <a href="#">Dissecting and Ligature Forceps</a>
<b>Aesculap</b>	BH139
<b>KLS Martin</b>	13-355-14

## AUTHORS

[DANDY](#)

## Medicon

16.11.84

## SPECIFICATIONS

[curved](#), [serrated](#)