

DANDY HAEMOSTATIC FORCEPS

curved; serrated

14.0cm (5.6")



SKU: 040-345-14

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Artery Forceps , Dissecting and Ligature Forceps |
| Aesculap | BH139 |
| KLS Martin | 13-355-14 |
| AUTHORS | DANDY |
| Medicon | 16.11.84 |
| SPECIFICATIONS | curved , serrated |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|--|
| UNIT | pcs. |
| TYPE | Artery Forceps , Dissecting and Ligature Forceps |
| Aesculap | BH139 |
| KLS Martin | 13-355-14 |

AUTHORS

[DANDY](#)

Medicon

16.11.84

SPECIFICATIONS

[curved](#), [serrated](#)