

DARDIK MULTI PURPOSE VASCULAR CLAMP

90° angled; atraumatic

13.5cm (5.3")




SKU: 142-022-90

Categories: [Open Surgery](#), [Thorax](#), [Cardiovascular Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vascular Clamps
Aesculap	FC314
KLS Martin	24-338-90
AUTHORS	DARDIK
SPECIFICATIONS	angled



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vascular Clamps
Aesculap	FC314
KLS Martin	24-338-90
AUTHORS	DARDIK

SPECIFICATIONS

[angled](#)