

## DARDIK MULTI PURPOSE VASCULAR CLAMP

30° angled; atraumatic

15.0cm (5.9")



**SKU:** 142-020-30

**Categories:** [Open Surgery](#), [Thorax](#), [Cardiovascular Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vascular Clamps</a>
<b>Aesculap</b>	FC302
<b>KLS Martin</b>	24-336-30
<b>AUTHORS</b>	<a href="#">DARDIK</a>
<b>SPECIFICATIONS</b>	<a href="#">angled</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vascular Clamps</a>
<b>Aesculap</b>	FC302
<b>KLS Martin</b>	24-336-30
<b>AUTHORS</b>	<a href="#">DARDIK</a>

## SPECIFICATIONS

[angled](#)