

## DARDIK MULTI PURPOSE VASCULAR CLAMP

15° angled; atraumatic

16.5cm (6.4")




**SKU:** 142-022-15

**Categories:** [Open Surgery](#), [Thorax](#), [Cardiovascular Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

|                       |                                 |
|-----------------------|---------------------------------|
| <b>UNIT</b>           | <a href="#">pcs.</a>            |
| <b>TYPE</b>           | <a href="#">Vascular Clamps</a> |
| <b>Aesculap</b>       | FC311                           |
| <b>KLS Martin</b>     | 24-338-15                       |
| <b>AUTHORS</b>        | <a href="#">DARDIK</a>          |
| <b>SPECIFICATIONS</b> | <a href="#">angled</a>          |



### ZUSÄTZLICHE INFORMATIONEN

|                   |                                 |
|-------------------|---------------------------------|
| <b>UNIT</b>       | <a href="#">pcs.</a>            |
| <b>TYPE</b>       | <a href="#">Vascular Clamps</a> |
| <b>Aesculap</b>   | FC311                           |
| <b>KLS Martin</b> | 24-338-15                       |
| <b>AUTHORS</b>    | <a href="#">DARDIK</a>          |

## SPECIFICATIONS

[angled](#)