

DARDIK MULTI PURPOSE VASCULAR CLAMP

45° angled; atraumatic

17.0cm (6.6")



SKU: 142-022-45

Categories: [Open Surgery](#), [Thorax](#), [Cardiovascular Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)


TYPE [Vascular Clamps](#)

Aesculap FC313

KLS Martin 24-338-45

AUTHORS [DARDIK](#)

SPECIFICATIONS [angled](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vascular Clamps](#)

Aesculap FC313

KLS Martin 24-338-45

AUTHORS [DARDIK](#)

SPECIFICATIONS

[angled](#)