

## DE BAKEY DIETHRICH BULLDOG CLAMP

straight; length: 50mm; jaw 14mm

closing force ca. 180g

**SKU:** 026-022-02

**Categories:** [Hemostats](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

|                       |  |
|-----------------------|--|
| <b>UNIT</b>           | <a href="#">pcs.</a>                                 |
| <b>TYPE</b>           | <a href="#">Bulldog Clamps</a>                       |
| <b>Aesculap</b>       | FB364  |
| <b>KLS Martin</b>     | 13-122-14  |
| <b>AUTHORS</b>        | <a href="#">DE BAKEY</a> , <a href="#">DIETHRICH</a> |
| <b>SPECIFICATIONS</b> | <a href="#">straight</a>                             |



### ZUSÄTZLICHE INFORMATIONEN

|                   |  |
|-------------------|--|
| <b>UNIT</b>       | <a href="#">pcs.</a>                                 |
| <b>TYPE</b>       | <a href="#">Bulldog Clamps</a>                       |
| <b>Aesculap</b>   | FB364  |
| <b>KLS Martin</b> | 13-122-14  |
| <b>AUTHORS</b>    | <a href="#">DE BAKEY</a> , <a href="#">DIETHRICH</a> |

## SPECIFICATIONS

[straight](#)