

DE BAKEY DIETHRICH BULLDOG CLAMP

angled; length: 44mm; jaw 10mm

closing force ca. 180g

SKU: 026-023-01

Categories: [Hemostats](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
Aesculap	FB363
KLS Martin	13-123-10
AUTHORS	DE BAKEY , DIETRICH
SPECIFICATIONS	angled



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
Aesculap	FB363
KLS Martin	13-123-10
AUTHORS	DE BAKEY , DIETRICH

SPECIFICATIONS

[angled](#)