

DE BAKEY DIETHRICH BULLDOG CLAMP

angled; length: 50mm; jaw 20mm

closing force ca. 180g


SKU: 026-023-03

Categories: [Hemostats](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
Aesculap	FB367
KLS Martin	13-123-20
AUTHORS	DE BAKEY , DIETRICH
SPECIFICATIONS	angled



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
Aesculap	FB367
KLS Martin	13-123-20
AUTHORS	DE BAKEY , DIETRICH

SPECIFICATIONS

[angled](#)