

DE BAKEY-SATINSKY BULLDOG CLAMP

angled; atraumatic

7.0cm (2.8")




SKU: 026-191-07

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
Aesculap	FB357
KLS Martin	24-291-07
AUTHORS	DE BAKEY , SATINSKY
SPECIFICATIONS	angled , atraumatic



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
Aesculap	FB357
KLS Martin	24-291-07
AUTHORS	DE BAKEY , SATINSKY

SPECIFICATIONS

[angled](#), [atraumatic](#)