

DE BAKEY SUCTION TUBE

removable suction tip; Luer connector

inner Ø 3.0mm; outer Ø 6.0mm; 27.0cm

SKU: 085-008-27

Categories: [Open Surgery](#), [Puncture / Special Needles](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Suction Tubes
Aesculap	GF946
AUTHORS	DE BAKEY
Medicon	09.85.07

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Suction Tubes
Aesculap	GF946
AUTHORS	DE BAKEY
Medicon	09.85.07