

## DERF-CONVERSE NEEDLE HOLDER

straight; longitudinal serrated

12.5cm (5")



**SKU:** 050-020-12

**Categories:** [Open Surgery](#), [Suture](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM200
<b>KLS Martin</b>	20-328-12
<b>AUTHORS</b>	<a href="#">CONVERSE</a> , <a href="#">DERF</a>
<b>Medicon</b>	10.12.82
<b>SPECIFICATIONS</b>	<a href="#">serrated</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>Aesculap</b>	BM200
<b>KLS Martin</b>	20-328-12

**AUTHORS**

[CONVERSE](#), [DERF](#)

**Medicon**

10.12.82

**SPECIFICATIONS**

[serrated](#), [straight](#)