

DIAMOND NITINOL LITHOTRIPTERBASKET

20 x 40mm; 8 wires

2.6mm; 230cm



SKU: PK-1LDXX2620230AK

Category: [Endoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Stonebaskets for Retrieval & Lithotripsy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Stonebaskets for Retrieval & Lithotripsy](#)