

DIETRICH BULLDOG CLAMP

straight; length: 50mm; jaw 12mm

closing force ca. 80g

SKU: 026-020-02

Categories: [Hemostats](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Bulldog Clamps](#)

Aesculap MB285

KLS Martin 13-120-10

AUTHORS [DIETRICH](#)

SPECIFICATIONS [straight](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Bulldog Clamps](#)

Aesculap MB285

KLS Martin 13-120-10

AUTHORS [DIETRICH](#)

SPECIFICATIONS

[straight](#)