

DIETRICH BULLDOG CLAMP

straight; length: 60mm; jaw 20mm

closing force ca. 80g

SKU: 026-020-03

Categories: [Hemostats](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--------------------------------|
| UNIT | pcs. |
| TYPE | Bulldog Clamps |
| KLS Martin | 13-120-12 |
| AUTHORS | DIETRICH |
| SPECIFICATIONS | straight |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--------------------------------|
| UNIT | pcs. |
| TYPE | Bulldog Clamps |
| KLS Martin | 13-120-12 |
| AUTHORS | DIETRICH |
| SPECIFICATIONS | straight |