

DIETHRICH BULLDOG CLAMP

angled; length: 40mm; jaw 8mm

closing force ca. 80g

SKU: 026-021-01

Categories: [Hemostats](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
KLS Martin	13-121-08
AUTHORS	DIETRICH
SPECIFICATIONS	angled



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
KLS Martin	13-121-08
AUTHORS	DIETRICH
SPECIFICATIONS	angled