

## DIETHRICH BULLDOG CLAMP

angled; length: 50mm; jaw 20mm

closing force ca. 80g


**SKU:** 026-021-03

**Categories:** [Hemostats](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bulldog Clamps</a>
<b>KLS Martin</b>	13-121-12
<b>AUTHORS</b>	<a href="#">DIETHRICH</a>
<b>SPECIFICATIONS</b>	<a href="#">angled</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bulldog Clamps</a>
<b>KLS Martin</b>	13-121-12
<b>AUTHORS</b>	<a href="#">DIETHRICH</a>
<b>SPECIFICATIONS</b>	<a href="#">angled</a>