

## DILATATOR 2 CHANNELS FOR INTRODUCTION

of guide wires for 21/24 Fr. sheath

300mm

**SKU:** 630-435-45

**Categories:** [Endoscopy](#), [Urology](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Nephroscopy</a>
<b>TYPE</b>	<a href="#">Dilatators</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Nephroscopy</a>
<b>TYPE</b>	<a href="#">Dilatators</a>