

## DILATATOR 2 CHANNELS FOR INTRODUCTION

of guide wires for 16.5/19.5 Fr. sheath

300mm

**SKU:** 630-435-40

**Categories:** [Endoscopy](#), [Urology](#)



### ZUSÄTZLICHE INFORMATIONEN

|              |                             |
|--------------|-----------------------------|
| <b>UNIT</b>  | <a href="#">pcs.</a>        |
| <b>RANGE</b> | <a href="#">Nephroscopy</a> |
| <b>TYPE</b>  | <a href="#">Dilatators</a>  |

### ZUSÄTZLICHE INFORMATIONEN

|              |                             |
|--------------|-----------------------------|
| <b>UNIT</b>  | <a href="#">pcs.</a>        |
| <b>RANGE</b> | <a href="#">Nephroscopy</a> |
| <b>TYPE</b>  | <a href="#">Dilatators</a>  |