

DOYEN INTESTINAL CLAMP

straight; diagonally serrated

21.0cm (8.4")



SKU: 160-210-21

Categories: [Open Surgery](#), [Stomach](#), [Intestines](#), [Rectum](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Stomach and Intestinal Clamp Forceps
Aesculap	EA124
KLS Martin	30-274-21
AUTHORS	DOYEN
SPECIFICATIONS	diagonally serrated , straight

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Stomach and Intestinal Clamp Forceps
Aesculap	EA124
KLS Martin	30-274-21
AUTHORS	DOYEN

SPECIFICATIONS

[diagonally serrated](#), [straight](#)