

## DOYEN MALLET

solid; head 35mm; 240g

21.0cm (8.4")


**SKU:** 191-212-35

**Categories:** [Open Surgery](#), [Orthopaedics](#)



### ZUSÄTZLICHE INFORMATIONEN

|                   |                         |
|-------------------|-------------------------|
| <b>UNIT</b>       | <a href="#">pcs.</a>    |
| <b>TYPE</b>       | <a href="#">Mallets</a> |
| <b>KLS Martin</b> | 23-204-03               |
| <b>AUTHORS</b>    | <a href="#">DOYEN</a>   |



### ZUSÄTZLICHE INFORMATIONEN

|                   |                         |
|-------------------|-------------------------|
| <b>UNIT</b>       | <a href="#">pcs.</a>    |
| <b>TYPE</b>       | <a href="#">Mallets</a> |
| <b>KLS Martin</b> | 23-204-03               |
| <b>AUTHORS</b>    | <a href="#">DOYEN</a>   |