

DOYEN SPHERICAL BURR

long Borchardt shaft

Ø 8.0mm



SKU: 190-290-08

Categories: [Open Surgery](#), [Orthopaedics](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Hand Drills
KLS Martin	22-210-08
AUTHORS	DOYEN

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Hand Drills
KLS Martin	22-210-08
AUTHORS	DOYEN