

## EAR NEEDLE

angled; 90° upwards; 0.3mm

16.0cm (6.4")

**SKU:** 111-648-03

**Categories:** [Open Surgery](#), [Otology](#)



## ZUSÄTZLICHE INFORMATIONEN

|                   |                         |
|-------------------|-------------------------|
| <b>UNIT</b>       | <a href="#">pcs.</a>    |
| <b>TYPE</b>       | <a href="#">Needles</a> |
| <b>KLS Martin</b> | 36-025-01               |



## ZUSÄTZLICHE INFORMATIONEN

|                   |                         |
|-------------------|-------------------------|
| <b>UNIT</b>       | <a href="#">pcs.</a>    |
| <b>TYPE</b>       | <a href="#">Needles</a> |
| <b>KLS Martin</b> | 36-025-01               |