

## ENDOOR C DISSECTING FORCEPS; INSERT

fenestrated; straight

Ø 5.0mm; 330mm



**SKU:** 611-210-53

**Categories:** [Endoscopy](#), [Laparoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

|              |                            |
|--------------|----------------------------|
| <b>UNIT</b>  | <a href="#">pcs.</a>       |
| <b>RANGE</b> | <a href="#">EndoorC</a>    |
| <b>TYPE</b>  | <a href="#">Dissectors</a> |

### ZUSÄTZLICHE INFORMATIONEN

|              |                            |
|--------------|----------------------------|
| <b>UNIT</b>  | <a href="#">pcs.</a>       |
| <b>RANGE</b> | <a href="#">EndoorC</a>    |
| <b>TYPE</b>  | <a href="#">Dissectors</a> |