

## ENDOOR C DISSECTOR; INSERT

straight; pitted; non-traumatic

Ø 5.0mm; 330mm



**SKU:** 611-200-53

**Categories:** [Endoscopy](#), [Laparoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">EndoorC</a>
<b>TYPE</b>	<a href="#">Dissectors</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">EndoorC</a>
<b>TYPE</b>	<a href="#">Dissectors</a>