

## ENDOOR C MIXTER GRASPING FORCEPS; INSERT

90° angled

Ø 5.0mm; 330mm

**SKU:** 611-205-53

**Categories:** [Endoscopy](#), [Laparoscopy](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">EndoorC</a>
<b>TYPE</b>	<a href="#">Dissectors</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">EndoorC</a>
<b>TYPE</b>	<a href="#">Dissectors</a>