

## ENGLISH MODEL TOWEL FORCEPS

sharp

9.0cm (3.6")

**SKU:** 030-109-09

**Categories:** [Hemostats](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Towel Forceps</a>
<b>Aesculap</b>	BF408
<b>KLS Martin</b>	13-877-09
<b>AUTHORS</b>	<a href="#">ENGLISCHES MODELL</a>
<b>Medicon</b>	17.50.28
<b>SPECIFICATIONS</b>	<a href="#">curved</a> , <a href="#">sharp</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Towel Forceps</a>
<b>Aesculap</b>	BF408
<b>KLS Martin</b>	13-877-09

## AUTHORS

[ENGLISCHES MODELL](#)

Medicon

17.50.28

## SPECIFICATIONS

[curved](#), [sharp](#)