

## ERCP-CATHETER; BOTTLE TOP SHAPED TIP

tube: 1.8mm/210cm

able for guide wire 0.035"

**SKU:** PK-3CXPk1810210M

**Category:** [Endoscopy](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [ERCP-Catheter, Papillotomes](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [ERCP-Catheter, Papillotomes](#)