

## ERCP-CATHETER; BOTTLE TOP SHAPED TIP

able for guide wire 0.035"

.8mm, length 250cm single use



**SKU:** PK-3CXEW1810250M-S

**Category:** [Endoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Instruments for Double-Balloon-Enteroscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Instruments for Double-Balloon-Enteroscopy](#)