

ERCP-CATHETER; SHORT-FILIFORM TIP

able for guide wire 0.035"

1.8mm, length 250cm single use



SKU: PK-3CTEW1810250-S

Category: [Endoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Instruments for Double-Balloon-Enteroscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

RANGE [Instruments for Double-Balloon-Enteroscopy](#)