

## ERCP-CATHETER; SHORT-FILIFORM TIP

able for guide wire 0.035"

1.8mm, length 250cm single use



**SKU:** PK-3CTEW1810250-S

**Category:** [Endoscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Instruments for Double-Balloon-Enteroscopy](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**RANGE** [Instruments for Double-Balloon-Enteroscopy](#)