

FERGUSSON (POPPEN) SUCTION TUBE

11 Fr.; 45° angled

14.5cm (5.8") working length



SKU: 085-034-11

Categories: [Open Surgery](#), [Puncture / Special Needles](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
TYPE [Suction Tubes](#)
AUTHORS [FERGUSSON, POPPEN](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
TYPE [Suction Tubes](#)
AUTHORS [FERGUSSON, POPPEN](#)